

Parameters & Features



Reusable / Rechargeable



Heart Rate



30-Day Wear



Phone Symptom Indicator



Water Resistant



No Cable Connection Required



Ultra Light and Flexible

Patient FAQs

Do I need to keep the mobile phone near me?

Yes, it is imperative to keep the provided phone within 10 feet of the patch and charge it at least every 12 hours to ensure continuous streaming of data.

Can I take a shower wearing my patch?

Yes, keep showers short with your back to the flow of water while showering. Do not use creams, shampoo, lotions, or similar products near the patch. After your shower, gently pat dry with a towel and minimize activity until the patch is dry.

What if my patch accidentally gets wet?

If the patch accidentally gets wet, gently pat dry with a towel and minimize activity until the patch is fully dry.

What do I do when monitoring is complete?

Remove the adhesive and discontinue wearing the monitor. Return ALL equipment to your Doctor's office or mail to Artella.



866-951-2142

713-821-3200

artellainc.com



Cardiac Device

Patient Information & Instructions



Your doctor has prescribed the **ARTELLA Solutions System** to monitor your heart. This is a wireless cardiac monitor that will be observed 24 hours a day. You will be able to maintain your normal daily activities since the patch is lightweight, discreet, and water-resistant.



Patient Device Instructions




- 1 To charge the patch, place it in the case, align with the gold dots, and wait for the light to turn green (fully charged).



- 2 Peel off strip 1 from the ECG disposable adhesive and align it with 3 dots on the back of the patch. Peel off adhesive strip 2 in order to place on the chest.




- 3 Common placement of the Arterra patch:
 - ▶ Measure 3 fingers below left clavicle.
 - ▶ Place at a 45 degree angle with point of heart facing the center of the chest.



Ensure skin is clean and free of hair with no lotions or oils applied where the device will be placed. Use alcohol swabs to clean the area and allow to dry before application.

Patient Phone Instructions

- 1 Turn on the phone by using the lower right side button. Tap on the Blue ARTELLA MCT icon.
- 2 Log in using the Monitoring Cycle and PIN code provided. The device battery percentage is displayed on the next page. The cloud indicates connection/streaming of your data.



- 3 If you have symptoms while wearing the device, tap on the picture on the screen that best represents what you are experiencing.

These devices are used as diagnostic tools only and are not life saving.
 For assistance, please call
866-951-2142 / 713-821-3200

To ensure accurate readings keep the following tips in mind:

- ▶ Make sure the provided phone remains on. Charge it nightly.
- ▶ Device battery is displayed on the phone app. Remove the ECG device from the chest and charge when the battery reaches 0%. Please charge during the DAY TIME.
- ▶ Stay within 10 feet of the phone. Keep the ARTELLA app installed.
- ▶ If the phone shows "No Patch" or "Poor signal", CALL ARTELLA.
- ▶ Avoid saunas and steam rooms and limit swimming. You can shower with the patch, but refrain from taking lengthy hot baths or showers.
- ▶ Do not use your nails or rub vigorously over the two outer circle sensors on the device.

When The Study is Complete:

- ▶ To remove the adhesive, slowly peel back one corner and blot the underlying area with water for a more comfortable removal. Do not pull on raised areas of the patch.
- ▶ Ensure the adhesive residue is completely off the device before reapplying a new adhesive by gently rubbing with the provided alcohol prep to remove it, on a flat surface.
- ▶ Remove and clean the patch, ensuring the residue is off, and place it in the white charging station. DO NOT THROW IT AWAY. Place the phone and charger in the ECG Monitoring Kit case that you received.
- ▶ Within 48 hours, return all items to your physician's office or by mail using the provided shipping label. FSA/HSA benefits may be accepted for any associated costs.